

SHOEBOX PureTest vs. SHOEBOX Standard/Pro: Key Differences

Feature	PureTest	SB Standard/Pro
Fleet Management - settings locked from Portal	✓ All test settings and workflow steps controlled from the web portal	✗ All test settings are set in the app
Enforce checks	✓ Functional check and Room Scan can be enforced to block testing	✗ Only the room scan can be enforced
Multi-user Login	✓ Each user may have a login under a single license	✗ Not available, only one user may be logged in at a time
Offline Shift Analysis	✓ Fully offline shift analysis available by selecting Projects for offline testing.	✗ Not available, must be online.
Workflow	✓ Wizard takes you through the entire testing process	✗ Steps must be selected to proceed through workflow
Single Results Document	✓ Results, baseline, shift analysis, hearing grade overlay and signatures can be viewed immediately after testing and does not require a data connection.	✗ A basic audiogram can be viewed immediately after testing. To view the summary report, it must be downloaded using a data connection.
Data Privacy	✓ Results are stored hidden until synced and then removed from the cache.	✓ Audiograms stay on the device until manually deleted. They are protected behind user logins and passcodes.
Languages	✗ Version 1 is English only	✓ Available in 27 subject-facing languages, and 4 Examiner-facing languages